



# ENGAUGE SAFETY Hypodermic Syringe with Safety Needle: Instructions For Use


The sterile TERUMO® ENGAUGE SAFETY SYRINGES for single use are intended for use in the aspiration and injection of fluids for medical purpose. After withdrawal of the needle from the body, the attached needle safety shield can be manually activated to cover the needle immediately after use to minimize risk of accidental needlesticks. The sterile ENGAUGE SAFETY Syringes for single use are single use, disposable devices, provided sterile.



**1.** Carefully remove the safety needle from the package. Attach the safety needle to a syringe with a push and clockwise twist motion.




**2.** Move the safety sheath away from the needle and toward the syringe barrel to the angle shown, prior to removing the protection cap. The fixing mechanism will keep the safety sheath in the position you set.




**3.** Remove the protection cap carefully to avoid damaging the needle. Perform injection procedure according to established techniques.


After the injection procedure, the safety mechanism could be activated by any of the three methods given below, i.e. 4a, 4b and 4c.




**STEP 4a.** Using a one-handed technique, push the tab forward with your thumb so that the sheath is less than 90 degrees from the needle. Then push down with a firm, quick motion until a distinct audible click is heard.



**STEP 4b.** Using a one-handed technique, push the tab forward with your finger so that the sheath is less than 90 degrees from the needle. Then push down with a firm, quick motion until a distinct audible click is heard.



**STEP 4c.** Using a one-handed technique, push the tab forward with your finger or thumb so that the sheath is less than 90 degrees from the needle. Position the sheath in the direction of activation at a lowest possible angle to the flat surface. Then push down with a firm, quick motion until a distinct audible click is heard.



**STEP 5.** Visually confirm that the needle is fully engaged under the lock. Dispose of used needles and materials following the policies and procedures of your facility as well as your federal and local regulations for sharps disposal.

**Note:** The location of the lock will vary depending upon needle size. Audible click may not be heard on small needle sizes; visual confirmation is required.

**RX ONLY.** The advertisement is directed to physicians only, and not to consumers. Refer to product labels and packaging insert for complete warnings, precautions, potential complications, and instructions for use.

## WARNINGS

- Pharmaceutical compounds prepared in syringes should be administered as soon as possible after they are removed from their original containers. As indicated in the ISO standard for syringes (ISO7886-1), syringe manufacturers are not required to perform testing for syringe compatibility with all possible injectable preparation.
- Visually inspect and carefully check the device and packaging before use. Improper transport and handling may cause structural and/or functional damage to device or packaging.
- Use once and discard immediately in accordance with local safety standards.
- Do not replace original needle cover after use.

## CAUTION

- Safety Needle is guaranteed non-toxic, sterile and non-pyrogenic if the individual package has not been opened or damaged.
- Do not attempt to deactivate the safety device by forcing the needle out of the safety sheath.
- For single use only. Do not reuse.
- Do not clean or re-sterilize.
- Safety Needle must not be reprocessed.
- If needle is bent or damaged, no attempt should be made to straighten needle or use the product.

## PRECAUTIONS

- Safety Needle should be used according to these instructions for use.
- Do not use if the individual package is damaged or opened.
- The product should be used immediately after opening of the package.
- Keep hands behind needle at all times during use and disposal.
- Store in a cool and dry place.
- Do not expose to heat or direct sunlight.
- Observe Universal Precautions on ALL patients.
- Do not use for high pressure injection of contrast media.

TERUMO and TERUMO are trademarks owned by Terumo Corporation, Tokyo, Japan, and they are registered with the U.S. Patent & Trademark Office. All rights reserved. All brand names are trademarks or registered trademarks of their respective owners. PM-07228