

# TERUMO® Conventional Needles & Syringes



## Superior sharpness and smoothness time after time

- Ultra Sharp Needles
- Clear Barrels
- Bold, Precise Scale Markings
- Smooth Plunger Action
- Tamper Evident Peel Packs


When comfort, performance and cost effectiveness are important, choose Terumo. Our ultra-thin wall needles are precisely honed and beveled for ease of penetration and less patient discomfort. Clear syringe barrels with bold scale markings permit quick visual inspection for precise dosage. An integrated plunger stop helps prevent the plunger from being accidentally pulled out. And our tamper evident blister packs result in minimal waste disposal.

### Terumo® Standard Hypodermic Needles

Packed: 100/box, 10 boxes/case (1,000)

Product Code	Description
■ NN1825R	18G x 1"
■ NN1838R	18G x 1½"
■ NN1925R	19G x 1"
■ NN1938R	19G x 1½"
■ NN2025R	20G x 1"
■ NN2038R	20G x 1½"
■ NN2125R	21G x 1"
■ NN2138R	21G x 1½"
■ NN2225R	22G x 1"
■ NN2238R	22G x 1½"
■ NN2325R	23G x 1"
■ NN2338R	23G x 1½"
■ NN2516R	25G x 5/8"
■ NN2525R	25G x 1"
■ NN2538R	25G x 1½"
■ NN2613R	26G x 1/2"
■ NN2713R	27G x 1/2"
■ NN2732R	27G x 1¼"
■ NN3013R	30G x 1/2"
■ NN3025R	30G x 1"

### All Competitors



- \* In Independent laboratory testing, Terumo needles were shown, on average, to be the sharpest among brands tested:  
**10% sharper than BD (22G - 30G)**  
**16% sharper than Covidien (22G - 27G)**

Call your Terumo representative at:

**1.800.888.3786**

or visit [www.TerumoTMP.com](http://www.TerumoTMP.com)

and click on Find Your Local Terumo Rep.

\* Data on file. Terumo Medical Products, April 2016.

These products have no components made of natural rubber latex.

When you think hypodermics...  
**Think sharp. Think Terumo.**



# Terumo® Conventional Needles & Syringes

## 1cc Tuberculin Syringes

**Graduations: 0.01cc (1cc total)**

TERUMO TB Syringes have clear, bold scale markings which assure dosage control and accuracy. The package is tamper evident, for greater protection of the sterile contents.

**Packed: 100/box, 10 boxes/case (1,000)**

### Product Code Description

SS-01T	1cc TB Syringe without needle, with luer slip
SS-01T2516	1cc TB Syringe with 25G x 5/8" removable needle
SS-01T2609	1cc TB Syringe with 26G x 3/8" removable needle
SS-01T2713	1cc TB Syringe with 27G x 1/2" removable needle



Catheter Tip



Luer Lock Tip



Eccentric Luer Tip



Luer Slip Tip

## 3cc Syringes with Needle

**Style: Luer Lock Tip**

**Graduations: 0.1cc (3cc total)**

**Packed: 100/box, 10 boxes/case (1,000)**

### Product Code Description

SS-03L2025	3cc, 20G x 1"
SS-03L2038	3cc, 20G x 1 1/2"
SS-03L2125	3cc, 21G x 1"
SS-03L2138	3cc, 21G x 1 1/2"
SS-03L2225	3cc, 22G x 1"
SS-03L2238	3cc, 22G x 1 1/2"
SS-03L2325	3cc, 23G x 1"
SS-03L2338	3cc, 23G x 1 1/2"
SS-03L2516	3cc, 25G x 5/8"
SS-03L2525	3cc, 25G x 1"

## 5cc Syringes with Needle

**Style: Luer Lock Tip**

**Graduations: 0.2cc (6cc total)**

**Packed: 100/box, 6 boxes/case (600)**

### Product Code Description

SS-05L2025	5cc, 20G x 1"
SS-05L2038	5cc, 20G x 1 1/2"
SS-05L2125	5cc, 21G x 1"
SS-05L2138	5cc, 21G x 1 1/2"
SS-05L2225	5cc, 22G x 1"
SS-05L2238	5cc, 22G x 1 1/2"

## 10cc Syringes with Needle

**Style: Luer Lock Tip**

**Graduations: 0.2cc (12cc total)**

**Packed: 100/box, 6 boxes/case (600)**

### Product Code Description

SS-10L2025	10cc, 20G x 1"
SS-10L2038	10cc, 20G x 1 1/2"
SS-10L2125	10cc, 21G x 1"
SS-10L2138	10cc, 21G x 1 1/2"

## 3cc Syringes without Needle

**Graduations: 0.1cc (3cc total)**

**Packed: 100/box, 10 boxes/case (1,000)**

### Product Code Description

SS-03L	3cc Luer Lock Tip
SS-03S	3cc Luer Slip Tip

## 5cc Syringes without Needle

**Graduations: 0.2cc (6cc total)**

**Packed: 100/box, 6 boxes/case (600)**

### Product Code Description

SS-05L	5cc Luer Lock Tip
SS-05S	5cc Luer Slip Tip

## 10cc Syringes without Needle

**Graduations: 0.2cc (12cc total)**

**Packed: 100/box, 6 boxes/case (600)**

### Product Code Description

SS-10L	10cc Luer Lock Tip
SS-10S	10cc Slip Tip
SS-10ES	10cc Eccentric Luer Tip

## 20cc Syringes without Needle

**Graduations: 1.0cc (25cc total)**

**Packed: 50/box, 4 boxes/case (200)**

### Product Code Description

SS-20L2	20cc Luer Lock Tip
SS-20S2	20cc Luer Slip Tip
SS-20ES2	20cc Eccentric Luer Tip

## 30cc Syringes without Needle

**Graduations: 1.0cc (35cc total)**

**Packed: 25/box, 8 boxes/case (200)**

### Product Code Description

SS-30L	30cc Luer Lock Tip
SS-30S	30cc Luer Slip Tip
SS-30ES	30cc Eccentric Luer Tip

## 60cc Syringes without Needle

**Graduations: 1.0cc (60cc total)**

**Packed: 25/box, 4 boxes/case (100)**

### Product Code Description

SS-60L	60cc Luer Lock Tip
SS-60ES	60cc Eccentric Luer Tip
SS-60C	60cc 2oz. Catheter Tip

## TerumoTMP.com

*These products have no components made of natural rubber latex.*

**RX ONLY.** The advertisement is directed to physicians only, and not to consumers. Refer to product labels and packaging insert for complete warnings, precautions, potential complications, and instructions for use. Products may not have regulatory approval in all countries. Please contact your local sales representative if you have questions about the availability of products in your area.

TERUMO and TERUMO are registered trademarks owned by Terumo Corporation, Tokyo, Japan. SurSaver is a registered trademark owned by Terumo Medical Corporation, Elkton, MD, U.S.A. PM-03444



TERUMO MEDICAL PRODUCTS