

Sur-Vet® Hypodermic Needles & Syringes



- Ultra Sharp Needles
- Clear Barrels
- Bold, Precise Scale Markings
- Smooth Plunger Action
- Tamper Evident Peel Packs

When comfort, performance and cost effectiveness are important, choose Terumo. Our thin wall needles are precisely honed and beveled for ease of penetration and less patient discomfort. Clear syringe barrels with bold scale markings permit quick visual inspection for precise dosage. An integrated plunger stop helps prevent the plunger from being accidentally pulled out. And our tamper evident blister packs result in minimal waste disposal.

Sur-Vet® Syringes with Needle

3cc Luer Slip, Packed: 100/box, 10 boxes/case (1,000)

Product Code	Description
■ 100090	20G x 1" RW
■ 100097	22G x 3/4" RW
■ 100098	22G x 1" RW
■ 100103	25G x 5/8" RW

3cc Luer Lock, Packed: 100/box, 10 boxes/case (1,000)

Product Code	Description
■ 100276	22G x 3/4" TW
■ 100277	22G x 1" TW
■ 100278	25G x 5/8" TW

1cc Luer Slip, Packed: 100/box, 10 boxes/case (1,000)

Product Code	Description
■ 100281	25G x 5/8" TW

Contact your Terumo representative at:
1.800.888.3786 or visit www.TerumoTMP.com
 and click on Find Your Local Representative


Sur-Vet® Hypodermic Needles

Packed: 100/box, 10 boxes/case (1,000)

Product Code	Description
■ 100204	20G x 1" RW
■ 100215	22G x 3/4" RW
■ 100279	22G x 3/4" TW
■ 100211	22G x 1" RW
■ 100220	25G x 5/8" RW
■ 100280	25G x 5/8" TW

RW - Regular Wall
 TW - Thin Wall

All Competitors



* In Independent laboratory testing, Terumo needles were shown, on average, to be the sharpest among brands tested:

10% sharper than BD (22G - 30G)

16% sharper than Covidien (22G - 27G)

Superior sharpness and smoothness time after time



* Data on file.

These products have no components made of natural rubber latex.

Sur-Vet® Hypodermic Needles & Syringes

Hypodermic Needles Packed: 100/box, 10 boxes/case (1,000)

Product Code	Description	Product Code	Description
■ NN1825R	18G x 1" TW	■ NN2338R	23G x 1½" TW
■ NN1838R	18G x 1½" TW	■ NN2525R	25G x 1" TW
■ NN1938R	19G x 1½" TW	■ NN2538R	25G x 1½" RW
■ NN1925R	19G x 1" TW	■ NN2613R	26G x ½" RW
■ NN2038R	20G x 1½" TW	■ NN2713R	27G x ½" RW
■ NN2125R	21G x 1" TW	■ NN2732R	27G x 1¼" RW
■ NN2138R	21G x 1½" TW	■ NN3013R	30G x ½" RW
■ NN2238R	22G x 1½" TW	■ NN3025R	30G x 1" RW
■ NN2325R	23G x 1" TW		



1cc Tuberculin Syringes

Graduations: 0.01cc (1cc total)

TERUMO TB Syringes have clear, bold scale markings which assure dosage control and accuracy. The package is tamper evident, for greater protection of the sterile contents.

Packed: 100/box, 10 boxes/case (1,000)

Product Code Description

SS-01T	1cc TB Syringe without needle, luer slip
SS-01T2516	1cc TB Syringe with 25G x 5/8" removable needle
SS-01T2609	1cc TB Syringe with 26G x 3/8" removable needle
SS-01T2713	1cc TB Syringe with 27G x 1/2" removable needle

3cc Syringes without Needle

Graduations: 0.1cc (3cc total)

Packed: 100/box, 10 boxes/case (1,000)

Product Code Description

SS-03L	3cc Luer Lock Tip
SS-03S	3cc Luer Slip Tip

20cc Syringes without Needle

Graduations: 1.0cc (25cc total)

Packed: 50/box, 4 boxes/case (200)

Product Code	Description
SS-20L2	20cc Luer Lock Tip
SS-20S2	20cc Luer Slip Tip
SS-20ES2	20cc Eccentric Luer Tip

5cc Syringes without Needle

Graduations: 0.2cc (6cc total)

Packed: 100/box, 6 boxes/case (600)

Product Code Description

SS-05L	5cc Luer Lock Tip
SS-05S	5cc Luer Slip Tip

30cc Syringes without Needle

Graduations: 1.0cc (35cc total)

Packed: 25/box, 8 boxes/case (200)

Product Code	Description
SS-30L	30cc Luer Lock Tip
SS-30S	30cc Luer Slip Tip
SS-30ES	30cc Eccentric Luer Tip

3cc Syringes with Needle

Style: Luer Lock Tip

Graduations: 0.1cc (3cc total)

Packed: 100/box, 10 boxes/case (1,000)

Product Code Description

SS-03L2025	3cc, 20G x 1"
SS-03L2038	3cc, 20G x 1½"
SS-03L2125	3cc, 21G x 1"
SS-03L2138	3cc, 21G x 1½"
SS-03L2238	3cc, 22G x 1½"
SS-03L2325	3cc, 23G x 1"
SS-03L2338	3cc, 23G x 1½"
SS-03L2525	3cc, 25G x 1"

10cc Syringes without Needle

Graduations: 0.2cc (12cc total)

Packed: 100/box, 6 boxes/case (600)

Product Code Description

SS-10L	10cc Luer Lock Tip
SS-10S	10cc Slip Tip
SS-10ES	10cc Eccentric Luer Tip

60cc Syringes without Needle

Graduations: 1.0cc (60cc total)

Packed: 25/box, 4 boxes/case (100)

Product Code Description

SS-60L	60cc Luer Lock Tip
SS-60ES	60cc Eccentric Luer Tip
SS-60C	60cc 2oz. Catheter Tip

RX ONLY. Refer to the product labels and package insert for complete warnings, precautions, potential complications, and instructions for use.

When you think hypodermics...
Think sharp. Think Terumo.



TERUMO and TERUMO are registered trademarks owned by Terumo Corporation, Tokyo, Japan. SUR-VET is a registered trademark owned by Terumo Medical Corporation, Elkton, MD, U.S.A. © 2022 Terumo Medical Corporation. All rights reserved. All brand names are trademarks or registered trademarks of their respective owners. PM-06431

