

SURFLO® Winged Infusion Sets

Terumo's soft, flexible interlocking wings readily conform to the body's contour. The wing is angled for ease of positioning of the needle during venipuncture, for less risk of vessel damage or occlusion of the lumen.

The precisely honed needle with thin wall design assures ease of penetration and continued access with a single insertion.

Product Code Description 3½" Tubing

Packed: 50/box, 10 boxes/case (500)

SV19BLSF	19G x ¾"	RW
SV21BLSF	21G x ¾"	TW
SV23BLSF	23G x ¾"	TW
SV25BLSF	25G x ¾"	RW

Product Code Description 8" Tubing

Packed: 50/box, 10 boxes/case (500)

SV25ELF	25G x ½"	RW
SV27ELF	27G x ½"	RW

Product Code Description 12" Tubing


Packed: 50/box, 10 boxes/case (500)

SV-18BLK	18G x ¾"	RW
SV-19BLK	19G x ¾"	RW
SV-21BLK	21G x ¾"	TW
SV-22BLK	22G x ¾"	TW
SV-23BLK	23G x ¾"	TW
SV-25BLK	25G x ¾"	RW

RW - Regular Wall

TW - Thin Wall

These products contain no components made of natural rubber latex.



Call your Terumo representative at:
1.800.888.3786 or visit www.TerumoTMP.com
and click on Find Your Local Terumo Representative

Compare our wings versus
the competition and see why
SURFLO® is a highly respected brand.

RX ONLY. The advertisement is directed to physicians only, and not to consumers. Refer to product labels and packaging insert for complete warnings, precautions, potential complications, and instructions for use. Products may not have regulatory approval in all countries. Please contact your local sales representative if you have questions about the availability of products in your area.

TERUMO and TERUMO are trademarks owned by Terumo Corporation, Tokyo, Japan, and they are registered with the U.S. Patent & Trademark Office. SURFLO is a registered trademark of Terumo Corporation, Tokyo, Japan and it is registered with the U.S. Patent & Trademark Office. © 2025 Terumo Medical Corporation. All rights reserved. All brand names are trademarks or registered trademarks of their respective owners. PM-06006V2